

## Sani-Cloth® 70

### Microbiological Efficacy Data: -

	Test Method	Organism	Contact Time	Result	Laboratory
Bacteria	EN16615	Pseudomonas aeruginosa	60 sec (c) & 2 min (d)	Pass, ≥Log 5	MGS
		Staphylococcus aureus	60 sec (c) & 2 min (d)		
		Enterococcus hirae	60 sec (c)		
	EN13727	Pseudomonas aeruginosa	30 sec (c)	Pass, ≥Log 5	Viroxy
		Staphylococcus aureus	30 sec (c)		
		Enterococcus hirae	30 sec (c)		
Yeast	EN16615	Candida albicans	60 sec (c) & 2 min (d)	Pass, ≥Log 4	MGS
	EN13624	Candida albicans	30 sec (c)	Pass, ≥Log 4	TecoLab
Enveloped Virus	EN14476	Vaccinia Virus	15 sec (c) & 60 sec (d)	Pass, ≥Log 4	Dr Brill & Partner (c) and Viroxy (d)
		Corona, HIV, Hep B & C, SARS, Flu	15 sec (c) & 60 sec (d)		
TB	EN14348	Mycobacterium terrae	60 sec (c)	Pass, ≥Log 4	Viroxy

(c) = Clean conditions

(d) = Dirty conditions

*\*Full Test reports available upon request*

Sarah Ashmore

Research & Development Manager, PDI

